



HOW TO CONFIGURE YOUR SYSTEM

CoolLED pE-1 excitation system

3 wavelength LED system with fixed dichroic mirror



3 Specify microscope. CoolLED provides adaptors to all current models.

2 CoolLED will provide a dichroic mirror at: 490nm or 515nm depending on wavelengths selected.

4 For additional LAMs and accessories www.CoolLED.com/Products/

Product code	Description	Quantity
244-1950	pE-1 Main Unit (Drives 3 λ)	1
244-1500	Control Pod	1
244-1600	pE-1 Fixed Collimator	1
244-1904	Light Guide (1m)	2
244-2200	Microscope Adaptor (Specify)	1
244-3400- λ_1 - λ_2	LAM @ $\lambda_1 + \lambda_2$ nm	1
244-3400- λ_3	LAM @ λ_3 nm	1

a **CoolLED** LED product
www.CoolLED.com