



HOW TO CONFIGURE YOUR SYSTEM

CoolLED pE-2 excitation system

4 wavelength LED system with exchangeable dichroic mirror



3 Specify microscope. CoolLED provides adaptors to all current models.

4 For additional LAMs and accessories www.CoolLED.com/Products/

| Product code | Description | Quantity |
|-------------------------------------|--------------------------------------|----------|
| 244-3750 | pE-2 Main Unit (Drives 4 λ) | 1 |
| 244-1500 | Control Pod | 1 |
| 244-4300 | pE-2 Exchangeable Collimator | 1 |
| 244-4312 | pE-2 Exchangeable dichroic mirror | 1 |
| 244-1904 | Light Guide (1m) | 2 |
| 244-2200 | Microscope Adaptor (Specify) | 1 |
| 244-3400- λ_1 - λ_2 | LAM @ λ_1 + λ_2 nm | 1 |
| 244-3400- λ_3 - λ_4 | LAM @ λ_3 + λ_4 nm | 1 |

a **CoolLED** LED product
www.CoolLED.com